

SENsYS_{EVO} TG-DSC

Simultaneous TGA and 3-D DSC



From -120 °C to 830 °C

- TG-DSC
- Coupling to gas analysis (BET, gas sorption, FTIR, MS)
- CALISTO software



The **SENSYS evo TG-DSC** offers the most precise DSC sensor together with a truly symmetrical microbalance to give **YOU** a combined TG-DSC system with totally unmatched performance. Additionally both measurement sensors are totally independent, so can be operated independently or simultaneously without compromising performance.

The HIGHLIGHTS

Differential Scanning Calorimeter


- **Incomparable precision:** the unique 3D DSC sensor totally surrounds your sample and reference crucibles so that the entire energy of any transformation is monitored thereby giving you an unequalled accuracy of measurement.
- **Parameter Independent Calibration:** because of the 3D transducer the calibration of **SENSYS evo** is independent of:
 - mass and form (powder, fibre, liquid, etc.) of the sample,
 - contact between the sample and the transducer,
 - crucible type,
 - sweeping gas (inert, oxidizing, reducing, wet, pressure) and flow rate.

Thermogravimetry

- **Excellent TG signal baseline:** the symmetry of the microbalance, with its sample and reference suspended vertically from the beam guarantee reproducible positioning of the sample over time and excellent stability of the TG signal baseline.
- **Automatic compensation of Buoyancy:** the sample and the reference are heated inside the same calorimetric block (symmetrical furnace), which means that the buoyancy effect is automatically compensated.

Ease of operation: dedicated market leading **CALISTO** software that is not only intuitive but powerful enough to perform every typical experiments and data treatment.

3D DSC SENSOR




3D Sensor

The sample and reference sensors are composed of **120 thermocouples mounted** in a cylinder that totally surrounds the measurement zone. These two cylinders can measure up to 94 % of all heat exchanged with the sample/reference, as compared to 20 % typically (50 % absolute maximum) with 2D place sensors.

The sensors are mounted in a calorimetric block that is further water cooled to eliminate any environmental variations and therefore giving you a highly precise and robust sensor with a unique level of specific sensitivity.

For studies at sub-ambient temperature, an automated cooling device utilizing evaporation of liquid nitrogen can be fitted onto the calorimetric unit.



Simultaneous TG & DSC

The performances of the **SENSYS evo TG-DSC** result from the **absolute symmetry** of this unique system in which the TG and DSC transducers, which are mechanically independent, retain their specific characteristics:

- **symmetry of the furnace:** the sample and the reference are heated inside the same calorimetric block, which means that the **buoyancy effect is compensated automatically**.
- **symmetry of the upper beam microbalance:** its vertical design with beams of identical length guarantees **perfectly** reproducible positioning of the sample over time and **excellent stability of the TG signal baseline**.


The **SENSYS evo TG-DSC** is an extremely sensitive instrument that can measure **actual variations in weight of the order of one microgram**.

Gas control panel

The **SENSYS evo TG-DSC** incorporates an automated gas control panel.

It is possible to:

- select from 3 different carrier gases (flow rate: 4 to 200 ml/min)
- mix one of these carrier gases with another "auxiliary or reactive" gaseous fluid (flow rate: 0.3 to 16 ml/min).



CRUCIBLES

Crucibles of 150 μ l capacity are available in aluminium, alumina, platinum or silica to suit the chemical nature of the sample. A silica tube can be added to protect the transducer and the analysis chambers during test involving sweeping with corrosive gas.



Evolved Gas Analysis

The **SENSYS evo** design features both a small furnace volume and also an easily accessible exhaust to allow for direct connection to further analysis devices such as MS, FTIR, GC etc. without excessive sample dilution or long transfer lines. Setaram offers a range of transfer lines and software solution for system comprehensive integration.

To date SENSYS has also been successfully coupled with:

- FTIR/MS and GC,
- Gas Sorption devices,
- An X-ray absorption spectroscopy apparatus.

Measurements under controlled humidity

Humidity is known to have a significant influence on the stability of many materials. In order to evaluate this specific influence, SETARAM has developed WETSYS a controlled humidity generator. It can:

- deliver a very precise level of relative humidity (RH) from 5 to 95 % depending on the temperature and measure the humidity around the sample.
- be connected to SETARAM instruments, and in particular to the **SENSYS evo TG-DSC** as well as other analysis instrument (RX, GC, etc.) for studies of sorption, hydration, stability, etc. See *SENSYS evo DSC brochure*



Applications

With its high precision and wide temperature range (-120 to 830 °C) the **SENSYS evo TG-DSC** has a wide range of applications, especially in the fields of thermodynamics (Specific Heat Capacity Cp precision within +/-1 %), pharmaceuticals (phase diagram, polymorphism, purity, thermal stability), safety of chemical process (transition, decomposition under high pressure), energy (catalysis, hydrogen storage, hydrogen adsorption), polymers (glass transition, measurement under pressure), etc.

View the application notes in your field, available for download, by visiting www.setaram.com!

A huge database is in the [Application Library](#) area of our website. We have also included a powerful search engine that will enable you to find the most applicable data.

Specifications

Temperature range	Ambient to 830 °C
With cooling accessories	-120 °C to 200 °C (Liquid Nitrogen Accessory)
Temperature accuracy	+/-0.1 °C
Temperature precision	+/-0.02 °C
Programmable temperature scanning rate (heating and cooling)	0.01 to 30 °C min ⁻¹
Cooling time	17 min (from 50 °C to -100 °C) with Liquid Nitrogen Accessory
DSC	
Enthalpy accuracy	+/-1 %
Calorimetric Precision	+/-0.1 %
RMS Noise	0.2 µW
Sensitivity (Joule effect at 30°C)	6.7 µV/mW
Resolution	0.35 µW / 0.035 µW
Dynamic range Resolution	+/-2500 mW
Gases	3 carrier gases (MFC from 4 to 200 ml/min) + 1 auxiliary or reactive gas (MFC from 0.3 to 16 ml/min)*
Crucibles	120 µl, 160 µl, 320 µl aluminium, incoloy, gold plated incoloy, graphite, alumina, platinum, etc.
Pressure (non controlled)	High Pressure Crucible (up to 500 bars / 7250 psi at 600 °C)
Pressure (measured & controlled)	High Pressure Crucible (up to 400 bars / 5800 psi at 600 °C)
TC	
Maximum balance capacity	3 g
Weight range	+/-200 mg
Weighing precision	+/-0.01 %
TG Resolution	0.02 µg, 0.002 µg
Baseline dynamic drift	< 10 µg
Isothermal temperature accuracy	+/-0.1 °C
Linear heating rate	30 °C/min
Vacuum	< 10 ⁻¹ mbar
Evolved gas	MS, FTIR, GC couplings
Weight	55 kg (121 lbs)
Dimensions (Height / Width / Depth)	60 closed, 80 opened / 53 / 58 cm (23,6 closed, 31.4 opened / 20.9 / 22.8 in)
Power requirements	230 V - 50/60 Hz

*MFC = Mass Flow Controller

Option: AKTS Thermokinetics software for comprehensive investigation of reaction or decomposition 



Offices in United Kingdom, Germany, Italy, Switzerland, China and Singapore

www.setaram.com
sales@setaram.com

SETARAM INSTRUMENTATION

7 rue de l'Oratoire
69300 Caluire - France
Phone +33(0)4 72 10 25 25
Fax +33(0)4 78 28 63 55

SETARAM Inc. (NJ office)

7905 Browning Rd, Suite 206
Pennsauken NJ 08 109-4323 - USA
Phone +1 (856) 910 9990
Fax +1 (856) 910 1150

SETARAM Inc. (HQ)

8430 Central Ave. Suite 3D
Newark, CA 94560 - USA
Phone +1 (510) 793 3345
Fax +1 (510) 402 4705