

## Analysis of Polyoxymethylene by TG-DSC + FTIR

### Introduction

One interesting feature of **LABSYS evo STA** is the possibility to connect its gas exhaust to a FTIR spectrometer in order to identify the vapours evolving from the sample and to monitor the level of each identified vapour.

### Experimental method

The following experimental procedure has to be used:


**Sample :** Polyoxymethylene

**Crucible :** Platinum

**Mass:** 41.37mg

**Atmosphere :** Argon 60ml/min

**Experimental procedure :** The temperature is programmed from ambient up to 600°C at 10 °C/min.



### Results and conclusions

The Heat Flow curve presents two endothermic effects. The low temperature peak (with a maximum at 187.2°C) may correspond to melting because it isn't related to any mass loss signal. While the higher temperature peak (with two maxima at 355.3°C and 378.9°C) corresponds to the decomposition of Polyoxymethylene and is related to the main mass loss of 98.82%. Finally, the last mass loss of 1.18% is due to a second decomposition.

From the FTIR spectra it is possible to draw the variation of intensity associated to the specific absorption band for carbon monoxide (2000-2200 $\text{cm}^{-1}$ ), formaldehyde (2600-3100 $\text{cm}^{-1}$ ) which are decomposition products of POM and tetrafluoroethylene (1100-1400 $\text{cm}^{-1}$ ). The release of this last compound is probably not linked with the degradation of the macromolecular chain, but with the presence of an additive or a residue of the preparation of the polymer in the final material.

**LABSYS evo**  
20°C to 1150°C  
20°C to 1600°C



[www.setaram.com](http://www.setaram.com)  
[sales@setaram.com](mailto:sales@setaram.com)

